



Community 1st
Care Network

Affordable Solutions for At Home Care Monitoring

Choices for independent living





Commissioning choices.....

- Need for efficiency savings within current system
- Opportunity for extended Primary Care services
- Appropriate care pathways
 - Re-design of referral patterns
 - Individual to patient requirements
- Efficient and accessible service provision





Community 1st
Care Network

The Community 1st Care Network

- Community 1st Care Network is a supply chain partner of Community 1st Oldham and is procurement-cleared to work with the local health community in Oldham
- We provide unique choices for reliable, easy to use and quick to install telecare and telehealth assistive technology
- These choices are solutions that support the delivery, or development, of tailored social and clinical care pathways

Choices for independent living





Telehealth and Telecare Benefits

- Reduction in hospital admissions
- Enables early discharge from hospital stay
- Efficient management of patients in the community
- Co-ordinated AHP monitoring environment
- Promotion of self-management of condition by client





Telecare and Telehealth Packages

- Enables patients to live comfortably and safely within their own home
- Offers an alternative to an extended hospital stay – monitor safely at home
- Supports living independently with 24/7 support if required
- Peace of mind for family and carers





Community 1st
Care Network

A Solution?



Close ByTM
NETWORK

Choices for independent living





Community 1st
Care Network

A Comprehensive Monitoring Solution



Choices for independent living





Community 1st
Care Network

Supporting Care at Home and Independent Living

CloseBy Network keeps those who need to know (the Care Network) notified and informed about your patients' activities of daily living



Choices for independent living



MEDICINE CABINET SENSOR



BED SENSOR



CHAIR SENSOR



FRONT DOOR SENSOR



KETTLE SENSOR




FRIDGE SENSOR



COOKER SENSOR



COMMUNITY 1ST CARE NETWORK'S SECURE SERVER




Choices for independent living





Community 1st
Care Network

Supporting Independent Living and Care at Home



Choices for independent living





Community 1st
Care Network

Real -Time Patient Monitoring (Activity Daily Living)



Community 1st
Care Network

Logged in as Yyyy | Log out

Monitor existing patient care pathway

Eric Wimp C1CN 001



Choices for independent living






Patient Assessment Tool



Logged in as Yyyyy | Log out



- Can the patient start their day without any help?
- Are they going to the bathroom too often?
- Are they sleeping in the chair during the day?
- Can they move about their home easily?






Community 1st
Care Network

Know What's Happening in a Patient's Home



Logged in as Yyyyyy | Log out




- Activities in the Past 24 Hours
- Status of Passive Sensors
- Quick View of Most Recent Activities

Choices for independent living





Know What's Happening in a Patient's Home




- Current Activities
- Past 24 Hours' Activities
- 7-Day Moving Average
- 30-Day Moving Average






Know What's Happening in a Patient's Home



- Activity Levels – View a 7-Day Moving Average




- or... over 30 Days





Know What's Happening in a Patient's Home




- Review individual sensor activity levels





Know What's Happening in a Patient's Home



- Graphic view of sequence of activities and events





Integrating Telehealth – Weight Scales

The screenshot displays the CloseBy Care Network interface. The top navigation bar includes 'Home', 'Sensor Status', and 'Sequence of Events'. The main content area is divided into three sections: 'Current Status', 'Recent Activity', and 'Total Activity'.

Current Status (highlighted with a red box):

Sensor Name	Status	Value
John Patient Weight	Online	187
John Patient Height	Offline	Offline
John Patient Fat Free Mass	Offline	Offline
John Patient Fat Ratio	Offline	Offline
John Patient Fat Mass Weight	Offline	Offline
Unassigned Data Weight	Offline	Offline
Unassigned Data Height	Offline	Offline
Unassigned Data Fat Free Mass	Offline	Offline
Unassigned Data Fat Ratio	Offline	Offline
Unassigned Data Fat Mass Weight	Offline	Offline

Recent Activity:

Time	Event
Today at 04:29:18 PM	Training Room Motion
Today at 04:29:10 PM	Training Room Motion
Today at 04:29:03 PM	Training Room Motion
Today at 04:28:59 PM	Training Room Motion
Today at 04:28:59 PM	Main Lobby Motion
Today at 04:28:59 PM	South Entry Door
Today at 04:28:47 PM	South Entry Door
Today at 04:28:43 PM	Main Lobby Motion
Today at 04:28:35 PM	Training Room Motion
Today at 04:28:29 PM	Training Room Motion
Today at 04:28:28 PM	Training Room Motion
Today at 04:28:28 PM	North Entry Door
Today at 04:28:27 PM	Humidity (Office)
Today at 04:28:24 PM	North Entry Door
Today at 04:28:22 PM	Training Room Motion

Total Activity (7 Day Average):

Day	Today (Events)	Yesterday (Events)
6	0	0
7	0	0
8	0	0
9	10	0
10	15	0
11	20	0
12	18	0
1	12	10
2	10	15
3	22	18
4	236	15
5	10	12
6	8	0
7	0	0
8	0	0
9	0	0
10	0	0
11	0	0
12	0	0
1	0	0
2	0	0
3	0	0
4	0	0
5	0	0






Activity Daily Living – Reports (24 Hours Summary)

CloseBy NETWORK		24-Hour Summary Report			24 Hour Summary Report - from 6:00am 05/15/11 to 6:00am 05/16/11																								
		Normal Within 30% of normal activity	Medium Change 31% - 60% of normal activity	Significant Change 61%+ of normal activity																									
ACTIVITY LEVEL	SENSOR NAME	ACTIVITY GRAPH			# OF EVENTS AVERAGE	MORNING HOURS				DAYTIME HOURS				NIGHTTIME HOURS															
		MORNING	DAY	NIGHT		6	7	8	9	10	11	12	1	2	3	4	5	6	7	8	9	10	11	12	1	2	3	4	5
●	Motion Sensor				<u>228</u> 237	10	18	26	4	6	10	2	4	10	22	6	13	33	16	26	10	8	15	9	3	2	1	4	2
●	Motion Sensor 2				<u>330</u> 270	3	23	52	10	24	18	2	24	16	40	10	4	11	28	33	20	8	15	7	3	1	4	1	
●	Front Door				<u>4</u> 3			2										1	1										
●	Kitchen Door				<u>14</u> 24			2	2					4	2				2	2									
●	Garage Upper Entry				<u>102</u> 67	1	5	18	2	10	10	2	4	4	12		2	8	14	6	4	3	2	1	1				
●	Garage Lower Entry				<u>10</u> 3			6											2									2	
○	Motion Sensor 2 Light Level				<u>34</u> NA	22	73	66	62	57	55	56	52	52	48	51	45	42	39	39	27	24							1





Activity Daily Living – Reports (Weekly Summary)





Community 1st
Care Network

Supported with our Unique Telehealth Device



Accurate remote monitoring of vital statistics

Choices for independent living





Community 1st
Care Network

Personal Medical File

Online

Medic4all

Safely preserves and stores personal medical history and examinations through telemonitoring

Allows consultation by physicians anywhere in the world (automatically translated in 18 languages)

Automatically feed data with examinations performed in diagnostic centers and private clinics

24/7 access to medical data via SMS-MMS or Internet for direct discussions with physicians



Choices for independent living





Value of Medical File



Health parameters monitoring results sent by Medic4all's monitoring devices

WristClinic™

Data received using the Wellbeing stations onsite



Relationship with affiliate Clinics:

- Referral
- Prescriptions
- Gyms
- Dentists



Physicians & specialists



Direct update of data with relevant personal and medical info by the final user





Community 1st
Care Network

Personal Medical File

Online

Medic4all

Medic4all

Language: English | Search: | Advanced Search

Welcome, Dr. Hershko | Logout

Home Patient Health Track Services My Account Mailbox Help

Home > Patient Health Track > Documents Manager

Documents Manager

Blood Pressure

Measurement Data	Time
150/90	25.11.2011 15:02
120/80	25.11.2011 15:02
150/90	25.11.2011 15:02
120/80	25.11.2011 15:02
150/90	25.11.2011 15:02
120/80	25.11.2011 15:02
150/90	25.11.2011 15:02
120/80	25.11.2011 15:02
150/90	25.11.2011 15:02
120/80	25.11.2011 15:02
150/90	25.11.2011 15:02
120/80	25.11.2011 15:02
150/90	25.11.2011 15:02
120/80	25.11.2011 15:02
150/90	25.11.2011 15:02
120/80	25.11.2011 15:02
150/90	25.11.2011 15:02
120/80	25.11.2011 15:02
150/90	25.11.2011 15:02
120/80	25.11.2011 15:02

Document Manager

- Measurements
- Respiratory Rate
- Blood Pressure
- Temperature
- Weight
- SPO2
- Glucose
- ECG2L
- ECG 1
- Consultations
- Medical History
- Medication Remider
- Nutrition & Diet

Last 10 Measurements

Show Measurements For: Any Time


Start Date: | End Date: |

Choices for independent living





Community 1st
Care Network



Choices for independent living





Community 1st
Care Network

Medic4all Doctor

Medic4all 

Health Maintenance Network

Home Portal / Marketplace

Language: English

Welcome tom111

Home HealthTrack My Account Mailbox (1) Patient Healthtrack Registration

tom111 PatientSearch

The current page has been customized from its template. [Revert to template.](#)

[Browse](#)

Patient Search

Patient Consultations

User Name:

First Name:

Last Name:

Date of birth:

NHS Number:

Gender:

Address:

Zip Code:



Choices for independent living





Community 1st

Care Network

Medic4all

How We Deliver

Platforms & Products



Choices for independent living





Self Monitoring Devices

NEW WristClinic™



- Monitor for blood pressure monitoring
- Heart rate and Heart rhythm regularity
- One Lead ECG
- Respiratory frequency
- Body temperature
- Sensor for Blood Oxygen Saturation (SpO2) monitoring
- Optional: accessory for constant monitoring of (SpO2)
- Possibility to integrate other medical devices through the serial port (RS-232)

New Wrist clinic has SMS and GPS inbuilt and touch screen

Open connectivity to other platforms for Providers or NHS

Particularlry for chronic and disease management at home

Can link with Medic4all Web Doctor or partners platform





Community 1st
Care Network

Contact us...

For further information or to arrange a demonstration contact our team.

www.community1stcarenetwork.co.uk

Community 1st Care Network Limited

Building 1000, Kings Reach, Yew Street, Stockport SK4 2HG

T: 0161 831 9722

E: info@community1stcarenetwork.co.uk

Choices for independent living

